



### Excellent sealing, Smooth inserttion

Unique Four-valve plus ultra-thin sealing membrane design, to ensure tightness and reduce the frictional resistance of instruments passing through



### More stable operation

The reinforced threaded tube design enhances the fixation of the abdominal wall, and the instruments do not slip or shift in and out



### Better visualization

Visible puncture rod, convenient to see the entire process of puncturing abdominal wall. Avoiding damage to the abdominal tissue

### Ordering information

Code	Diameter	working length
AVTC 5 - 100	5	100
AVTC 10 - 100	10	100
AVTC 12 - 100	12	100
AVTC 15 - 100	15	100

- Sinolinks Medical Innovation, Inc.
- Add.: No.10 Ziwei Road, Zhonglou Zone, Changzhou, Jiangsu, China, 213023
- Tel:0519-88021690
- E-mail:sales@sinolinks.com.cn
- Web site:www.sinolinks.com.cn

# Disposable Trocar



### Application

Used for minimally invasive surgical procedures to establish a path of entry for endoscopic instruments for primary and secondary insertions.